

COMPATIBLE ITEMS

# GloTran™ Hydrogen Peroxide Plasma Disinfectant



A Better Way to  
Destroy Pathogens

GloTran™ was designed to reduce Hospital Acquired Infections caused by Noninvasive, Portable Clinical Items and Non-Critical Equipment in healthcare. Manual cleaning takes a significant contact time with the device, while also exposing staff to potentially hazardous chemicals. For items that are difficult to clean effectively in a manual process, GloTran™ can get the job done quickly and safely. GloTran's advantages over other methods of disinfection include low temperature, dry treatment environment and fast cycle times of 20 minutes or less.



### Items compatible for disinfection in GloTran™\*

#### COMMONLY USED NON-CRITICAL MEDICAL EQUIPMENT

Airflow Sensors (sleep lab)  
 Bag-Valve-Mask (BVM) Devices  
 Bladder scanner – Handheld  
 Blood coagulation meter  
 Blood Pressure Cuffs  
 Bone Conduction Oscillator  
 Breast pump  
 Clippers/ tweezers  
 CPAP hoses  
 Defibrillator paddles  
 Dermatoscopes  
 Ear Gauges  
 Ear Impression Syringes  
 Ear Lights (used for otoblock insertion)  
 ECG cables and lead wires  
 Electric Thermometers  
 Fingernail care equipment  
 Gait belts  
 Headbands  
 Headphones  
 Hearing aids  
 Hearing aid stethoscope  
 Humidifier parts  
 Infusion pumps  
 Immittance emission probe tips  
 Limb Props  
 Mammography paddles  
 Measuring tape  
 Nebulizer mouthpieces or facemasks  
 Needle exchange tongs  
 Non-critical suction equipment  
 Ophthalmoscopes  
 Orthopedic Tourniquets  
 Otoacoustic emission probe tips

Oto-lights  
 Otosopes + Otoscope specula  
 Oximeter probe  
 Peak flow meters (Respiratory equipment)  
 Probe tubes for real-ear measurement  
 Prosthetics  
 Pulse Oximeters  
 Reflex Hammers  
 Roho cushions  
 Safety glasses or goggles  
 Scissors (for non-surgical purposes)  
 Sequential compression device sleeves  
 Speculum lights  
 Spirometer  
 Splints  
 Stethoscope + Case  
 Swallowing assessment kits  
 Tuning forks  
 Ultrasound Transducers

#### COMMONLY USED CLINICAL ITEMS

Call Buttons  
 Cell Phones, pagers,  
 Computer keyboards and mice  
 I.D. Badges  
 Pens  
 Remotes  
 Tablets (iPads, etc.)

... and hundreds of other non-critical care and clinical tools regularly used in healthcare facilities.

### Items NOT currently compatible with GloTran™

#### SURGICAL TOOLS, LIFE-SUPPORTING AND CRITICAL CARE DEVICES

Airway support devices  
 Dissecting pins  
 Endoscopes  
 Gowns or linens  
 Laryngoscopes and blades  
 Scalpels  
 Scissors  
 Surgical Forceps

Other surgical equipment requiring >99.999% sterility

#### OR ITEMS COMPOSED OF THE FOLLOWING MATERIALS

Wood, Cotton, Paper, Brass, Zinc, Iron, Silver, Bronze, Chrome, Copper, Chromium, Manganese, Liquid or powder.

\*NOTE: Items which are grossly soiled or contain visible contamination must be cleaned before use in GloTran.