

FASTEST/MOST EFFECTIVE WAY TO COMBAT SKIN INFECTIONS!

PROVIDE THE SAFEST POSSIBLE ENVIRONMENT FOR COMPETITION

THE STERILASER DELIVERS 57,000 MICROWATTS OF GERMICIDAL ULTRAVIOLET UVC LIGHT, DISTROYING 99.5% OF ALL HARMFUL BACTERIA, FUNGI, AND VIRUSES IN LESS THAN 0.25 SECONDS. ALL COMMON PATHOGENS FOUND ON MATS THAT CAUSE DISEASES SUCH AS RINGWORM, IMPETIGO, AND MRSA ARE KILLED IN ONE PASS OF THE STERILASER. WHEN PULLED AT A NORMAL WALKING PACE, IT COMPLETELY STERILIZES AN ENTIRE 40' X 40' WRESTLING MAT IN 4 MINUTES



KEY FEATURES

4 POWERFUL GERMICIDAL LAMPS



47" WIDTH/OPTIONAL BATTERY PACK



SHATTERPROOF COATED



Compare for Yourself

Time to Disinfect

Entire Mat

Frequency of Use

Knowledge Level

Effectiveness



STERILASER

4 Minutes

Can be used between matches or rounds - as frequently as needed or

desired.

Can be used by anyone walking at a normal pace

Eliminates 99.7% of pathogens with

one pass

When you're done

Easily rolls into storage closet

MOP & BUCKET

45 Minutes

Can only be applied once before a match or tournament.

Proper use includes knowledge of the solution, application and disposal. insufficient knowledge of the process makes the application useless.

Virtually ineffective if not applied correctly.

Dumping mat cleaner down the drain is in violation of federal EPA law.

Proper disposal must be done in an EPA approved manner.

CONTACT INFO

Phone: 877-904-7111 Email: highlandsfightgear@gmail.com

